Composito Universal hibrido fotopolimerizavel reforçado com Nanopartículas
Composite universal hibrido fotopolimerizable reforzado con Nanopartículas
Nanoparticle-reinforced photopolymerizable hybrid Universal Composite

PORTUGUÊS

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ENGLISH

1. Identification of the Product

FILL Magic is a hybrid composite made of fine particles, with a high load and reinforced with nanoparticles. FILL Magic is radiopaque and photopolymerizable in the dental cavity or in molds.

1.1. Indications

- Class I, II, III, IV and V direct reconstructions
- Aesthetic corrections
- Semidirect restorations including inlays, onlays, and facets

1.2. Technical data

Average particle size: 0.5 µm. Average filler percentage by weight: approximately 75%. Area of application: universal (anterior and posterior). Handling: excellent consistency, high condensability, highly shapeable and excellent color stability.

1.3. Composition

Bis-GMA, Bis-EMA, UDMA, TEGMA, photoinitiator, fillers and pigments.

1.4. Presentation

Replacement

1 syringe containing 4g in the colors:

Enamel: A1, A2, A3, A3.5, A4, B1, B2, B3, C1, C2, C3, C4, D3

Dentin: A1, A2, A3, B2 and C2

Incisal

01 Instruction for use

2. Application Techniques

- Carry out prophylaxis of the tooth to be treated.
- Select the appropriate shade of the material for restoration.
- Implement absolute isolation.
- Prepare the cavity according to the recommendations for composite restorations.
- Implement pulp protection, when necessary, with calcium hydroxide or glass ionomer type liner material according to the clinical assessment.
- Used the acid conditioning technique according to the manufacturer's instructions.

- Apply the adhesive according to the manufacturer's instructions.
- Place the necessary amount of FILL Magic inside the cavity preparation with a suitable instrument, adjusting it firmly to the cavity walls.
- Light-cure each layer for 40 seconds. After removal of the matrix the composite should be subjected to additional polymerization.
- Work around the surface of the restoration with stones or fine diamond tips for finishing.
- Shape the proximal surface with abrasive strips for finishing and polishing.
- Check the occlusion with articulating paper. Carefully adjust the occlusion by removing the material with fine-polished diamond tip or multilaminate tip.
- Polish with mounted silicone tips or disks and abrasive strips for polishing and finishing.

3. Precautions, Restrictions and Warnings

3.1. Precautions

- This product is manufactured to be used only as stated in the instructions for use. Usage other than described in the instructions for use is the sole and exclusive responsibility of the professional.
- This resin is light sensitive. Therefore, excessively bright environments can cause premature polymerization.
- Always keep the product closed when not in use.
- After being ejected outside of the packing it must be used quickly.
- Keep the product away from children.
- It should not be used with products containing (phenol, clove oil, copal varnish) or zinc oxide cements with eugenol, as these may prevent or affect the polymerization.

- Discoloration may occur when using cationic mouthwashes, as well as plaque and chlorhexidine developers.
- After use the material must be properly closed.
- All dental material must be used solely by dentists and dental laboratories.

3.2. Restrictions and Warnings

There is no contraindication when FILL Magic is used by a dental professional in accordance with the Instructions for Use. However, it is not recommended to be used in patients who have a history of allergy to the components of this product.

4. Validity

The expiration date shown on the packaging represents the maximum period for use of the product. Under proper storage conditions, FILL Magic is valid for 36 months. Do not use after the expiry date.

5. Storage

It must remain in its original packaging and be closed after use. Store at room temperature (15°C - 30°C) in a dry place, protected from light and heat.

Direct light on the product may change the properties of FILL Magic.

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