# Magic BOND DE com Flúor

## PORTUGUÊS

### Descrição

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## ENGLISH

### Description

Magic Bond DE with fluorine is a bonding agent light-curing adhesive with filler and fluoride for enamel and dentin (Primer and Bond).

Magic Bond DE is compatible with all restorative materials: primers, light curing composite resin and resin cement.

### Composition

Urethane dimethacrylate, pyromellitic dimethacrylate (PMDM), 2-hydroxyethyl methacrylate, bisphenol A glycidyl methacrylate, Anhydrous, 1,3-glycerol dimethacrylate, 2,4,6-trimethylbenzoyldiphenyl phosphine oxide, camphorquinone, butylated hydroxytoluene, Silicon Dioxide, bis (hydroxyethyl)-aminopropyl N-hydroxyethyl octadecylamine dihydrofluoride.

### Indications

Direct restorations of composites, light curing comonomers and comonomers. Adhesive cementation of indirect restorations: ceramic, ceromer and composite (veneers, inlays, onlays and crowns).

### Application Technique

#### Direct Restorations

1. Apply Magic Acid - a 37% phosphoric acid gel - on the enamel and dentin. The phosphoric acid should stay in contact with the enamel for at least 15 seconds, and in contact with the dentin during, maximum, 10 seconds.

2. Rinse the acid with water.

3. Air dry enamel and dentin, avoiding dentin dehydration, once the dentin surface should stay slightly humid.

4. Apply Magic Bond DE with fluorine over the enamel and dentin, by brushing it gently for 10 seconds and, thereafter, apply a stream of air to disperse the material and get a thin layer. Avoid very thin layers.

5. Light-cure Magic Bond DE with fluorine for 20 seconds. The surface should look glossy

and shiny, otherwise, repeat the application. 6. There after, apply a direct restorator material, such as Fill Magic, Fill Magic NT Premium, Concept Advanced, Brilliant or other similar material according with the indicated instructions make the finish restoration.

### Indirect Restoration

1. Prepare the restoration surface of ceramic, ceromer or composite, according to the product's use instructions. Restorations should be prepared by making an acid conditioning of the surfaces according to the material and to the required technique.

2. Apply Magic Bond DE with fluorine all over the surfaces of the restoration and, spreading the material with a stream of air, until a very thin and uniform layer is reached. At this stage, do not light cure the adhesive, keeping it out of light.

3. Apply Magic Acid - a 37% phosphoric acid gel - on the enamel and dentin. The phosphoric acid should stay in contact with the enamel for at least 15 seconds, and in contact with the dentin during, maximum, 10 seconds.

4. Rinse the acid with water.

5. Air dry enamel and dentin, avoiding dentin dehydration, once the dentin surface should stay slightly humid.

6. Apply Magic Bond DE with fluorine over the enamel and dentin, by brushing it gently for 10 seconds and, thereafter, apply a stream of air to disperse the material and get a thin layer. Avoid very thin layers.

7. Light-cure Magic Bond DE with fluorine for 20 seconds. The surface should look glossy and shiny, otherwise, repeat the application.

8. Cement the restoration by using Fill Magic Dual Cement, according to the instructions.

### Contraindications

There are no contraindications when Magic Bond DE with fluorine is used by a dentist

according to the manufacturer's instructions. Nevertheless, in those patients who has allergy to any of its components, we advise not to use it, or yet when the indicated technique cannot be applied.

### Precautions

All dental material should only be used by dentists and authorized laboratory technician. Keep the product away from reach of children.

Close the bottle cap immediately after use. Do not store together with material containing eugenol or oil of cloves, which may inhibit Magic Bond DE with fluorine polymerization.

### Validity and Storage

Magic Bond DE with fluorine is suitable for its specific end application, if stored at environment temperature (15°C - 30°C), and to the shelter of sunlight. The validity date which is in the package represents the maximum period the product should be used.

### Presentation

01 bottle of Magic Bond DE with fluorine with 5 mL.

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**Bonding agent with light curing charge and one-component for enamel and dentin**