EDTA TRISSÓDICO LÍQUIDO

PORTUGUÊS

Descrição

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ENGLISH

Description

LIQUID EDTA 18% is a chelating agent used inside the pulp cavity to prepare the walls of root canals. It acts to demineralize dentin particles, facilitating their dissolution and efficiently removing dentin sludge and scrapings, increasing dentin permeability and facilitating adherence of the treated connective tissue to the root surface in periodontal therapy.

Presentation:

01 bottle with 20mL

01 Instructions for use

Composition:

Disodium Ethylenediaminetetraacetic Acid, Sodium Hydroxide and Deionized Water.

Indication:

- Facilitate endodontic instrumentation of root canals;

- Removal of the smear layer (dentin mud);

- Preparing the dentin wall for the adhesion of obturating materials, as well as improving the quality of core cementation by removing smear from the inside of prepared canals.

- It acts as a demineralizer of dentin particles by chelating calcium and magnesium ions.

Instructions for use:

In endodontics.:

1. irrigate the root canal with 1 to 2 drops of 18% LIQUID EDTA;

2.Leave the product to act for 2 to 3 minutes to begin instrumentation.

3.The reaction with the dentin takes 5 to 6 minutes, after which time the product becomes less effective and neutralizes itself;

4.Repeat the canal irrigation procedure, applying the same amount of product as necessary.

5.Keep the canal moistened with 18% LIQUID EDTA throughout the procedure;

6.At the end of the instrumentation, irrigate the canal using an irrigating solution such as sodium hypochlorite.

In prosthetics:

1. Once the canal has been prepared, place a few drops of 18% LIQUID EDTA inside the canal and leave it to act for 5 minutes;

2. Irrigate with an irrigating solution such as sodium hypochlorite;

3. Repeat if necessary;

4. Dry the canal with absorbent paper tips and cement according to the technique indicated by the manufacturer of the chosen product.

Precautions:

- Completely remove 18% LIQUID EDTA before root canal restoration or temporary dressing

- Isolate the area before using the product.

- In case of direct contact with oral mucosa, rinse with plenty of water.

- In case of contact with eyes, rinse thoroughly with water. Seek medical help for evaluation and subsequent treatment.

Warning:

- FOR PROFESSIONAL USE

- Respect the expiration date. Never use the product if it has expired.

- Keep out of reach of children.

- Empty packaging must be disposed of in accordance with the rules for disposing of other materials, products and packaging from dental practices/laboratories.

Contraindications

- There are no contraindications when LIQUID EDTA 18% is used by a dental professional in accordance with the instructions for use. However, in patients who have a history of allergy to any of the components, we recommend not using it, or when the stipulated technique cannot be applied.

Shelf life:

Under proper storage conditions, LIQUID EDTA 18% has a shelf life of 36 months, even after opening, provided it is kept in its original packaging, properly closed and packaged according to the storage instructions.

The expiry date on the packaging represents the maximum period of use for the product. Do not use after the expiration date.

Close the bottle immediately after use to avoid oxidation of the 18% LIQUID EDTA.

Storage:

- Keep the bottle tightly closed immediately after use.

- Protect from moisture, heat and direct sunlight.

- Store the product at temperatures between 15°C and 30°C.

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