# PASTA ENDO PTC

## PORTUGUÊS

**Descrição**

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## ENGLISH

**Description**  
ENDO PTC is a water-soluble lubricant gel based on carbamide peroxide which, together with sodium hypochlorite, releases oxygen, facilitating the loosening of pulp fragments and the penetration of cleaning instruments, helping to remove dentin residues from the pulp cavity during chemical-mechanical root canal preparation and preventing darkening of the pulp tissue. Facilitates the sliding of endodontic files and prevents them from breaking inside the canals.

**Presentation:**  
ENDO PTC - REPLACEMENT  
01 syringe with 1g; 2g; 3g; 4g; or 5g.  
02 applicator tips  
01 instructions for use  
ENDO PTC - KIT  
02; 03; 04; or 05 syringes, with contents of 1g; 2g; 3g; 4g; or 5g.

02 applicator tips for each syringe  
01 Instruction for use - printed or digital.  
**Composition:**  
Carbamide peroxide, Polysorbate 80, Polyethylene glycol.  
**Indication:**  
Aid in the chemical-mechanical preparation of the root canal, to be used alternately with Sodium Hypochlorite solutions (Dakin's Liquid or Milton's Liquid).

**Technical information:**  
Each component of ENDO PTC has its own function:  
Polyethylene glycol is the carrier and serves as a transport for the substances, making their association viable. The vehicle of the formulation is a greasy substance with a detergent capacity that is completely soluble in water and stable even at room temperature. This characteristic of being a greasy substance enables it to act as a lubricant during instrumentation, allowing it to run smoothly and effectively, maintaining the necessary pace. It thus helps to reduce the risk of instruments getting stuck in the canal walls, preventing fractures or other accidents.  
Polysorbate 80 is a detergent that acts by reducing surface tension, allowing the product to spread more quickly and evenly over the entire surface of the root

canal, favoring the cleaning process by facilitating the penetration of substances into the dentinal tubules. It also has moisturizing and emulsifying properties, adsorbing grease and keeping it in suspension, making it easier to remove. Carbamide peroxide reacts with sodium hypochlorite, resulting in a reaction between them. This reaction promotes the release of a large amount of oxygen, in a process of effervescence which, by mechanical action, removes the remains of dentin excised by the instruments, organic substances, bacteria and toxic products from the root canal to the outside environment.

**Instructions for use:**  
1.After performing the pulpectomy and determining the working length, irrigate the pulp chamber abundantly with Sodium Hypochlorite (e.g. Dakin's or Milton's liquid).  
2.Fill the pulp chamber with ENDO PTC, select the appropriate file (manual or rotary), carry out instrumentation to remove dentin debris.  
3.Irrigate abundantly with sodium hypochlorite (e.g. Dakin's or Milton's liquid).  
4.Apply ENDO PTC as necessary.  
5.Rinse the pulp chamber thoroughly with water to completely remove the residue and irrigation solution.

6.In cases of retreatment, ENDO PTC can be used in conjunction with the solvent to help remove the end of the gutta-percha cone, following the same root canal preparation procedures.

**Precautions:**  
- Completely remove ENDO PTC before root canal restoration or temporary dressing  
- Isolate the area before using the product.  
- In case of direct contact with oral mucosa, rinse with plenty of water.  
- In case of contact with eyes, rinse thoroughly with water. Seek medical help for evaluation and subsequent treatment.

**Warning:**  
- FOR PROFESSIONAL USE  
- Respect the expiration date. Never use the product if it has expired.  
- Keep out of reach of children.  
- Empty packaging must be disposed of in accordance with the rules for disposing of other materials, products and packaging from dental practices/laboratories.components are advised not to use it, or when the stipulated technique cannot be applied.

**Contraindications:**  
- There are no contraindications when ENDO PTC is used by a dental professional in accordance with the instructions for use. However, patients with a history of allergy to any of the components are advised not to use it, or when the stipulated technique cannot be applied.

**Shelf life:**  
Under proper storage conditions, ENDO PTC has a shelf life of 24 months, even after opening, provided it is kept in its original packaging, properly closed and packaged according to the storage instructions.

The expiry date on the packaging represents the maximum period of use for the product. Do not use after the expiration date.  
Close the syringe immediately after use to avoid oxidation of ENDO PTC.

**Storage:**  
- Keep the syringe tightly closed immediately after use.  
- Protect from moisture, heat and direct sunlight.  
- Store the product at temperatures between 5°C and 30°C.

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