# FILL Magic FLOW

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## ENGLISH

### 1. Product Description

Fill Magic Flow is a low viscosity microhybrid composite that can be applied to prepared dental cavities without the need for pressing or packing allowing a strong bond to the cavity prepared walls.

After curing, its properties are similar to those of a regular hybrid composite. It can also be applied in small quantities to characterize pits, mask stains and discolorations, which could influence the aesthetic result. Can also be used as an opaquer.

### 1.1. Indications

- Dental enamel restoration with erosion, abrasion or movement damage.
- Application of aesthetic layer of enamel surface to hide defects and restores.
- Sealing of pits and fissures.
- Filling of areas with difficult access.
- Class V.
- Small restorations of classes I, III and IV.
- Direct facet resin.
- Repairs in porcelain.
- Cementation of porcelain veneers.
- Dental fragments collage.
- Restoration of structural defects: hypoplastic spots.
- Marginal defects repair in composite resin restorations.
- Construction of small nuclei.
- Characterization and opacification of the tooth structure.

### 1.2. Composition

Bis-EMA, TEGMA, BHT, EDAB, Camphoroquinone, fillers and pigments.

### 1.3. Presentation

Fill Magic Flow Refill: 01 syringe with 1,2g in shades: Opaque; A1; A2; A3; A3,5; B2; C3 / 02 application cannulas / 01 Instruction for use

### 2. Method of Application

Prepare the cavity according to the usual technique. Use an adhesive such as Magic Bond, according to its usage instructions. Apply Fill Magic Flow directly inside the cavity, using the application cannula, in thin layers (2mm maximum), for optimizing the light-curing.

After the use, pull the syringe plunger a few millimeters, to avoid product overflow. Light cure for 40 seconds. Use the conventional methods and materials for sculpting. Polish with non-abrasive disks or tips.

### 3. Precautions, Restrictions and Warnings

#### 3.1. Precautions

- This product is manufactured to be used only as set forth in the Instructions for use. The use differently than described in the operating instructions are the sole and exclusive responsibility of the professional to use.
- This resin is light sensitive. So too intense ambient light can cause premature polymerization.
- Keep the product always closed when not in use. After being ejected out of the package, it should be used quickly.
- Keep out of the reach of children.
- It should not be used with products containing phenol, clove oil, copal varnish or zinc oxide eugenol cements, as they prevent or affect the polymerization.
- Discoloration can occur when used cationic mouthwashes and plaque revealing and chlorhexidine.
- After using the material, it must be properly closed.
- Every dental material should be used only by dentists and dental laboratories authorized.

#### 3.2. Restrictions and Warnings

There are no restrictions when Fill Magic Flow is used by dental professionals according to the usage instructions. However, it is not recommended for patients with a history of allergy to any of its components.

#### 4. Expiry

The expiry date printed on the packaging represents the maximum date for using this product. Do not use after the expiry date. Under adequate storing conditions, Fill Magic Flow has a shelf life of 36 months.

#### 5. Storing

Must be kept in its original packaging and closed after each use. Keep at room temperature (15°C - 30°C), in a dry place, protected from light and heat. Direct light over the product may alter the properties of Fill Magic Flow.

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