PORTUGUÊS

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ENGLISH

1. Product identification
Adsil is an addition-type silicone-based impression material for use in dentistry. It comes in Putty Soft (high consistency), Light Body (low viscosity) and Regular Body (medium viscosity) versions.
Adsil has excellent dimensional stability, high elastic recovery capacity and great reproduction of details.

1.1. Indications for use
1.2. Adsil Putty Soft
Material for the simultaneous or double working impression technique
Molding of dentate, partially dentate or edentulous patients, for the realization of prostheses on implants, fixed or removable, partial or total prostheses and veneers.
1.3. Adsil Light Body / Adsil Regular Body
Material for the simultaneous or double working impression technique
Molding of dentate, partially dentate or edentulous patients, for the realization of prostheses on implants, fixed or removable, partial or total prostheses and veneers
Used as a precision material in the second stage of impression taking

1.4. Composition
Base: Silicon dioxide, Organopolysiloxanes, Surfactant, Pigments and Flavors.
Catalyst: Silicon dioxide, Organopolysiloxanes, Pigments and Platinum agent.

2. Instructions for use
2.1 Clinical Times

	Putty Soft
Mixing time	0'30"
Working time	2'00"
Prey time	3'30"

Tipo	Light Body/Regular Body Bisnaga/Tube
Mixing time	0'30"
Working time	2'00"
Prey time	3'30"

Tipo	Light Body/Regular Body Cartridge
Mixing time	0'10"
Working time	1'00"
Prey time	3'30"

2.2 Double Molding
Adsil Putty Soft
1. Select a rigid tray capable of holding 3 to 4 mm thick of impression material;
2. Clean and dry the field and, if necessary, remove the gums using a gingival retractor wire;
3. Using a measuring spoon, remove equal quantities of Adsil Putty Soft base and Adsil Putty Soft Catalyst;
4. Use your fingertips to mix the base with the catalyst for 30 seconds to obtain a mixture with a uniform and homogeneous color; do not wear latex gloves.
5. Remember not to use latex gloves when handling the silicone;
6. Place the mixture in the tray and take the first impression, according to the time indicated in the clinical time table;
7. Before molding, make a relief with PVC film or any other plastic.
8. If you do not choose this type of relief, after inserting the tray into the patient's oral cavity, wait an average of 1 minute to provide relief through movements;
9. Remove the tray from the oral cavity and reserve the tray with the impression to receive the fluid base;
10. Relief can also be carried out after the impression has been taken, using a scalpel.

Cartidge - Adsil Light Body / Adsil Regular Body
1. insert the Adsil Light Body or Adsil Regular Body cartridge into the dispenser.
2. Remove the cartridge cap and fit the mixing tip
3. Pulse the dispenser until the product comes out;
4. Discard the first output of product, to ensure uniformity in the dispersed quantity of base and catalyst in the mixture;
Note: With each use of the product or change of the mixing tip, a quantity of product must be dispensed, thus guaranteeing an output of equal quantities of base and catalyst.
5. Take part of the material to the cavity preparation.

Note: It is advisable to adapt the intraoral tip to the mixing tip in order to apply the product directly to the cavity preparation. This technique can ensure a mold free of bubbles or porosity.
6. Remove the plastic film from the impression with Adsil Putty Soft, if the relief has been made.
7. If necessary, to ensure a better impression, disperse a small amount of the mixture inside the first impression, i.e. the reserved Adsil Putty Soft impression;
8. Position the tray containing the Adsil Putty Soft impression in the desired location, exerting pressure for 4 to 5 seconds and keeping it stable until it is removed from the oral cavity;
9. After removing the tray, rinse the mold with running water and disinfect before pouring the model;
10. After use, remove the mixing tip and return with the original cartridge cap. The mixing tip and/or intraoral tip must not be reused.

Tube - Adsil Light Body / Adsil Regular Body
1. On a glass plate or spatulation block, spread equal amounts of Adsil Light Body or Adsil Regular Body base and catalyst, enough to make the desired molding.
Note: Spreading a greater or lesser amount of Adsil Light Body Catalyst or Adsil Regular Body Catalyst directly affects working and setting times, i.e. a greater amount speeds up setting, while a lesser amount slows down setting;
2. Using a spatula, mix the base with the catalyst, spatulating vigorously for 30 seconds in order to homogenize all the material;
3. Transfer the mixture to a molding syringe and take part of the material to the cavity preparation;
4. Remove the plastic film from the impression with Adsil Putty Soft, if the relief has been done.
5. Disperse the rest of the material inside the first impression, i.e. the reserved Adsil Putty Soft impression;
6. Position the tray containing the Adsil Putty Soft impression in the desired location, exerting pressure for 4 to 5 seconds and keeping it stable until it is removed from the oral cavity;
7. After removing the tray, rinse the mold with running water for 30 seconds and disinfect before pouring the model.

2.3 Simultaneous impression taking
Adsil Putty Soft + Tube Adsil Light Body or Adsil Regular Body
1. Select a rigid tray capable of holding 3 to 4 mm thick impression material;
2. Clean and dry the field, if necessary, remove the gums using a gingival retractor;
3. With the aid of a measuring spoon, remove equal quantities of Adsil Putty Soft Base and Adsil Putty Soft Catalyst;
4. Simultaneously, on a glass plate or spatulation block, another professional and/or assistant should disperse equal quantities of Adsil Light Body or Adsil Regular Body Base and Catalyst, enough to make the desired impression.
Note: Spreading a greater or lesser amount of Adsil Light Body Catalyst or Adsil Regular Body Catalyst directly affects working and setting times, i.e. a greater amount speeds up setting, while a lesser amount slows down setting;
5. Using the fingertips, mix the Adsil Putty Base with the Adsil Putty Catalyst for 30 seconds, while at the same time the other professional and/or assistant, using a spatula, should start mixing the Base with the Catalyst, spatulating vigorously for 30 seconds, in order to homogenize all the material to obtain a mixture with uniform colouring; do not wear latex gloves.

6. Transfer the Adsil Putty Soft mixture to the tray, while at the same time another professional and/or assistant transfers the Adsil Light Body or Adsil Regular Body mixture to an impression syringe;
7. Apply the Adsil Light Body or Adsil Regular Body mixture contained in the impression syringe to the cavity preparation, immediately positioning the tray containing the Adsil Putty mixture in place, applying pressure for 4 to 5 seconds and keeping it stable until it is removed from the oral cavity;
8. After removing the tray, rinse the tray with running water for 30 seconds and disinfect before pouring the model.
Note: In this procedure, no relief is performed
Adsil Putty Soft + Cartucho - Adsil Light Body / Adsil Regular Body
1. Select a rigid tray with capacity to provide a thickness of 3 to 4 mm of impression material;
2. Clean and dry the field; if necessary, remove the gums using a gingival retractor;
3. Using a measuring spoon, take equal amounts of Adsil Putty Soft base and Adsil Putty Soft Catalyst;
4. Using your fingertips, mix the base with the catalyst for 30 seconds to obtain a uniformly colored mixture and transfer the mixture to the tray; Do not wear latex gloves.
5. At the same time, another professional and/or assistant should fit the Adsil Light Body or Adsil Regular Body cartridge into the dispenser, remove the cartridge cap and fit the mixing tip, pulsing the dispenser containing the Adsil Light Body or Adsil Regular Body cartridge until the product comes out through the tip;
6. Discard the first spout of product to ensure uniformity in the amount of base and catalyst dispersed in the mixture;
Note: Each time the product is used or the mixing tip is changed, a quantity of product must be dispensed to ensure that the base and catalyst come out in equal quantities.
7. Apply the mixture to the cavity preparation or under the mixture of Adsil Putty Soft with catalyst in the tray, immediately positioning the tray in place, applying pressure for 4 to 5 seconds and keeping it stable until it is removed from the oral cavity;
Note: It is advisable to adapt the intraoral tip to the mixing tip to apply the product directly to the cavity preparation. This technique can ensure a mold free of bubbles or porosity.
8. After removing the tray, rinse the tray with running water for 30 seconds and disinfect before pouring the model.
Note: No relieving is performed in this procedure.
3. Disinfection
- After removal from the mouth, rinse the impression under running water for 30 seconds, disinfecting it with a suitable commercially available dental disinfectant solution that does not affect the surface or dimensions of the impression.
- Acrylic trays must be protected against water absorption.
- If hydrogen peroxide is used for disinfection, it is advisable to rinse the impression with plenty of warm water to prevent bubbles from forming.
4. Cleaning the trays
- The material should not be removed using a sharp instrument.
- Residues can be removed by soaking in commercially available alcohol or acetone and rinsing under running water.
Only use these solvents in well-ventilated areas.
- Clean and disinfect the trays using the usual methods.
- If hydrogen peroxide is used for disinfection, it is advisable to rinse the impression with plenty of warm water to prevent bubbles from forming.

4. Cleaning the trays
- The material should not be removed using a sharp instrument.
- Residues can be removed by soaking in commercially available alcohol or acetone and rinsing under running water.
Only use these solvents in well-ventilated areas.
- Clean and disinfect the trays using the usual methods.
- If hydrogen peroxide is used for disinfection, it is advisable to rinse the impression with plenty of warm water to prevent bubbles from forming.

5. Precautions, warnings and contraindications.
5.1. Precautions
- Avoid direct exposure of the material to sunlight or other sources of intense light.
- Avoid changing the base and catalyst caps, as polymerization may occur.
- Avoid contact with skin and eyes.
- In the event of accidental skin contact, wash immediately with soap and water and seek medical advice in the event of allergies or local discomfort.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

5.2. Warning
- All dental materials should only be used by authorized dentists and dental laboratories.
- Keep out of reach of children.
- This product is intended to be used only as described in the instructions for use. Any use other than that described in the instructions for use is the sole and exclusive responsibility of the professional using it.
- Always wear gloves when working, except latex gloves when handling Adsil.
Skin secretions, latex gloves and latex-contaminated surfaces can interfere with the polymerization process of polyvinylsiloxanes.
- Materials and surfaces that the impression may come into contact with (countertops, sinks, teeth, preparations, retraction cords, etc.) must be free of latex residues.
- Wear vinyl/nitrile gloves, or if this is not possible, wash the gloves with warm water and detergent before handling the product and/or the impression.
- Compounds containing eugenol or hemostatic substances can also impede the perfect polymerization process.
- Do not pour the plaster until 1 hour after casting and, for the best quality model, use type IV plaster.

5.3. Contraindications
There are no contraindications when Adsil is used by a dental professional in accordance with the instructions for use. However, we recommend not using it on patients with a history of allergy to any of its components, or when the stipulated technique cannot be applied.

6. Shelf life
- Adsil is valid for 36 months from the date of manufacture, even after opening, provided the storage conditions are respected.
- Do not use after the expiry date.
- The expiry date on the packaging represents the maximum period of use of the product.

7. Storage
- Store at 15°C to 30°C, in a dry place, protected from light and heat.
- After use, close the packaging tightly and immediately.
- Keep the product in its original packaging.
- Improper storage of the product can lead to degradation or malfunction, directly affecting shelf life.

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