TRUE CORE

PORTUGUÊS

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ENGLISH

1.1 Description

TrueCore is a multipurpose, dual (photo and self-polymerizable), radiopaque resin cement, indicated for cementing bridges, inlays, onlays and composite crowns, ceramics (feldspar, lithium disilicate, leucite, alumina, zirconia, hybrid ceramics), alloys metal and metal-ceramic, in addition to fixed prostheses, adhesive prostheses, intra-radicular retainers (in metallic alloys or fiberglass) and creation of coronal cores or fillings associated or not with intra-radicular retainers (pins). The product is sold in a auto-mixing double-body syringe, which makes it easier to obtain the correct proportion of each paste (1:1), facilitating clinical work.

1.2 Composition

Base: Bis-GMA, Bis-EMA, TEGDMA, BHT, EDAB, camphorquinone, barium and aluminum borosilicate, silicon dioxide and pigments.

Catalyst: TEGDMA, Bis-EMA, Bis-GMA, BHT, benzoyl peroxide, silicon dioxide and barium aluminum borosilicate.

1.3 Indications

TrueCore is indicated for:

Cementation with the associated use of an adhesive system for bridges, inlays or onlays made of ceramics (feldspar, lithium disilicate, leucite, alumina, zirconia, hybrid ceramics), metallic alloys, metal-ceramics or composites and other hybrid materials.

Cementation with associated use of adhesive system for full crowns, fixed prostheses, adhesive prostheses, intra-radicular retainers.

Manufacture of coronal cores or fillings associated or not with intra-radicular retainers, works carried out using any of the materials used for these described prosthetic purposes.

1.3 Presentation

Refill of 6g

1 double-body syringe of TrueCore containing 6g, 3g of base + 3g of catalyst

Available in Yellow colors; Opaque; Universal or Bleach

4 mixing tips

4 applicator tips

1 instruction for use

Kit

4 double-body syringes of TrueCore, containing 24g, 3g of base + 3g of catalyst in each syringe.

This can be: 1 syringe of each color, or if of the same colors, depending on the colors available (Yellow; Opaque; Universal or Bleach)

16 mixing tips

16 applicator tips

1 instruction for use

2. Instruction for Use

2.1 Application Techniques

Treatment of prosthetic parts

Metallic or metal-ceramic parts:

- Clean the prosthetic part with a 70% alcohol solution to remove any manufacturing residues such as dust and grease.

- Dry the piece well with compressed air free of oil and moisture.

- Blast the inner surface of the part with aluminum oxide particles (50 µm) at moderate pressure (approximately 2-2.5 bar).

- This process creates micro-roughness on the surface, increasing the contact area and improving the mechanical retention of the resin cement.

- Apply an MDP-based metallic primer.

Ceramic pieces based on zirconia or alumina:

- Clean the alumina/zirconia prosthetic part with a 70% alcohol solution to remove any manufacturing residues, such as dust and grease.

- Dry the piece well with compressed air free of oil and moisture.

- Blast the inner surface of the part with aluminum oxide particles (50 µm) at moderate pressure (approximately 2-2.5 bar).

- This process creates micro-roughness on the surface, increasing the contact area and improving the mechanical retention of the resin cement.

- For zirconia ceramics: Apply an MDP-based zirconia primer.

- Apply a silane agent to the inner surface of the zirconia/alumina part. The silane acts as a chemical bonding agent between the ceramic and the resin cement.

- Let the silane dry for a few minutes (according to the manufacturer's instructions) and avoid touching the treated surface.

Ceramic parts such as Lithium Disilicate (E-max), Lithium Silicate with Zirconia, Lithium Disilicate with Aluminosilicate:

- Clean the prosthetic part with a 70% alcohol solution to remove any manufacturing residues such as dust and grease.

- Dry the piece well with compressed air free of oil and moisture.

- Etch with 5 to 10% hydrofluoric acid for 20 seconds

- Wash thoroughly with jets of air and water.

- Dry the piece well with compressed air free of oil and moisture.

- Apply a silane agent to the inner surface of the part. The silane acts as a chemical bonding agent between the ceramic and the resin cement.

- Let the silane dry for about 1 minute (or according to the manufacturer's instructions) and avoid touching the treated surface.

- After this process, a thin layer of adhesive can be applied, this step being optional.

Feldspathic, Leucite and Polymer Infiltrated Ceramic Network System (PICN) parts:

- Clean the prosthetic part with a 70% alcohol solution to remove any manufacturing residues such as dust and grease.

- Dry the piece well with compressed air free of oil and moisture.

- Etch with 10% hydrofluoric acid for 60 seconds

- Wash thoroughly with jets of air and water.

- Dry the piece well with compressed air free of oil and moisture.

- Apply a silane agent to the inner surface of the part. The silane acts as a chemical bonding agent between the ceramic and the resin cement.

- Let the silane dry for about 1 minute (or according to the manufacturer's instructions) and avoid touching the treated surface.

- After this process, a thin layer of adhesive can be applied, this step being optional.

Cementation of indirect restorations

- Remove the provisional.

- Completely eliminate the fixation cement from the cavity preparation and carry out prophylaxis.

- Isolate the operating field.

- Properly clean the preparations with a pumice stone and cup for prophylaxis, washing and drying thoroughly with a jet of air and water.

- Dry the dental element well with compressed air free of oil and moisture.

- Condition the enamel (for 30 seconds) and the dentin for (15 seconds), if a conventional adhesive is used, with 37% phosphoric acid.

- Wash with plenty of water.

- Apply the preparation with absorbent paper, leaving the dentin tissue slightly moistened.

- Apply the Magic Bond DE adhesive actively with the Technobrush micro-appliator and light-activate it

according to the manufacturer's instructions.

- In areas or situations where light from the light-curing device cannot reach, TrueCore cement is not compatible with acidic pH adhesives (simplified acid-etching adhesives and some single-step self-etching adhesives). The high acidity of these adhesives compromises the chemical curing of the resin cement at the interface with the adhesive. In these cases you can choose two possibilities:

- Use of a Dual adhesive system on the tooth surface

- Apply a thin layer of hydrophobic adhesive (Magic Bond) over the adhesive layer of Magic BOND DE and light-activate (It is important to ensure that you do not apply a thick layer that would harm the laying of the piece; to do this, remove the excess with a clean, dry brush the part during preparation and remove before light curing the hydrophobic adhesive). In both options, any incompatibility problem between adhesive and cement will be resolved.

- Select the appropriate cement color Double-body syringe: Connect the self-mixing tip to the double-body syringe, discard the first portion and fill the prosthetic part with a sufficient amount of cement.

- Apply the cement inside the prosthetic piece, place it over the preparation, pressing gently, remove excess with a brush and dental floss, or over the tooth cavity covering the entire surface.

- Wait approximately 5 minutes for chemical curing to begin. - After removing excess, light cure the TrueCore for approximately 40 to 60 seconds on each side.

Cementation of Intra-root Posts: (metallic or fiberglass):

- Perform prophylaxis on the dental element with a pumice stone, water and a prophylaxis cup, washing and drying thoroughly with a jet of air and water.

- Isolate the operative field and unblock the conduit following the biomechanical principles of intra-radicular retainers.

- Clean the inside of the canal with a brush to remove any traces of endodontic cement or filling material.

- Wash and dry well with oil- and moisture-free compressed air.

- Etch enamel (for 30 seconds) and dentin (for 15 seconds) with 37% phosphoric acid. The inside of the conduit must be conditioned with specific tips.

- Wash with plenty of water.

- Dry the canal well with absorbent paper cones, leaving the dentin tissue slightly moistened.

- Apply the Magic Bond DE adhesive actively with the Technobrush micro applicator and carry out the light activation according to the manufacturer's instructions.

- In areas or situations where light from the light-curing device cannot reach, TrueCore cement is not compatible with acidic pH adhesives (simplified acid-etching adhesives and some single-step self-etching adhesives). The high acidity of these adhesives compromises the chemical curing of the resin cement at the interface with the adhesive. In these cases you can choose two possibilities:

- Use of a Dual adhesive system on the tooth surface

- Apply a thin layer of hydrophobic adhesive (Magic Bond) over the adhesive layer of Magic BOND DE and light-activate (It is important to ensure that you do not apply a thick layer that will harm the laying of the piece; to do this, remove excess with a clean brush or tip absorbent paper and place the pin inside the conduit and remove it before light-activating the hydrophobic adhesive for at least 60 seconds).

- In both options, any incompatibility problem between adhesive and cement will be resolved.

- Select the appropriate color of cement.

- Connect the self-mixing tip to the double-body syringe, discard the first portion and fill the prosthetic part with a sufficient amount of cement.

</