PORTUGUÊS

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ENGLISH

1. Product Description

NT Premium Flow is a nanohybrid composite resin with high fluidity, designed to be applied to prepared dental cavities, without

the need for pressure or condensation, creating a close bond on the walls of the preparation. After curing, its properties resemble those of well-known hybrid composite resins. In small amounts it can also be applied to characterize grooves, pigments or to mask dental staining or discoloration that may affect the esthetic result, as well as opacifying restorative surfaces in composite resins or metal surface repairs. Indication Dental restoration of enamel with erosion, abrasion or dislocation damage. Application of an esthetic layer on the enamel surface in order to hide defects and restorations. Pit and fissure sealing. Filling in undercut areas. Class V. Small classes I, III and IV restorations. Direct resin facet. Repairs in porcelain. Cementation of porcelain veneers. Bonding of dental fragments. Restoration of structural defects: hypoplastic stains. Repairs to marginal defects in composite resin restorations. Construction of small cores. Characterization and Opacification of the dental structure.

2. Indication

- Due to its lower viscosity, it is more easily applicable in class III and V restorations, and is also indicated for minimally invasive cavity preparations (including small occlusal restorations and seals that are not subject to tension) and the repair of marginal defects (filtration) in previously existing restorations.

- It is also indicated for repairs to small esthetic defects in indirect restorations and provisional resin or acrylic ones.

- Owing to its fluidity and elasticity (lower level of fillers), it is ideal for the liner and base (initial layers) of class I and II restorations.

3. Presentation Reposition

- 01 syringe containing 2g in the colors: NT Premium WO; NT Premium TR; NT Premium A1; NT Premium A2; NT Premium A3; NT Premium A3,5; NT Premium A4; NT Premium B2; NT Premium B3; NT Premium C2.

- 02 Cannulas for application.

- 01 Instruction for use.

4. Composition

Monomers (Bis-GMA, Bis-EMA, TEGMA), BHT, photoinitiators, fillers and pigments.

5. Instructions for use

- Prepare the cavity according to the usual technique.
- Use a Magic Bond DE type adhesive as per the instructions for use.

- Apply NT Premium Flow directly with the dispensing tip, into the cavity.

- Do not apply layers more than 2 mm thick in order to optimize polymerization.

- Once the application is complete, pull back the syringe plunger to avoid extra material escaping.

- Light cure for 40 seconds.

- Use the conventional method and materials for sculpting.

- Use non-abrasive disks or tips for polishing.

6. Precautions, Restrictions and Warnings

6.1. Precautions

This product is manufactured to be used only as stated in the instructions for use. Usage other than described in the instructions for use is sole and exclusive responsibility of the

professional.

Contact with saliva and blood during composite application may damage the restoration. We recommend the use of an absolute isolation made from cotton or rubber protection or other suitable isolation technique.

This resin is light sensitive. Therefore, excessively bright environments can cause premature polymerization. Keep the product closed when not in use. After being injected outside of the packing it should be used quickly. Keep the product away from children. It should not be used with products containing (phenol, clove oil, copal varnish) or zinc oxide eugenol cements, as these may prevent or affect the polymerization.

Discoloration may occur when using cationic mouthwashes, as well as plaque and chlorhexidine developers. After use the material must be properly closed. All dental material must be used solely by authorized dentists and dental laboratories.

6.2. Restrictions and Warnings

There are no contraindications when NT Premium Flow is used by a dental professional according to the instructions for use.

However, it is not recommended to be used in patients who have a history of allergy to the components of this product. As a preventive measure to avoid a possible reaction in the pulp, calcium hydroxide cement, glass ionomer cement or another type of cavity liner may be used prior to insertion of material into the cavity.

7. Expiry Date

Under appropriate storage conditions, NT Premium Flow has a validity of 36 months. The expiration date shown on the packaging represents the maximum period of use of the product. Do not use after the expiry date.

7. Storage

It must remain in its original packaging and be closed after use. Store at room temperature (15° - 30°C), in a dry place, protected from light and heat. Direct light on the product can change the properties of NT Premium Flow. If the product is stored under refrigeration, it should be removed from the refrigerator before use in order to return to ambient temperature.

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